# Senso








Use it. Leave it. Fold it.

promedicare.eu



### **About**

The SEAT & BACK Positioning Systems of the SENSO Range switch from an unfolded-use configuration to a folded non-use configuration avoiding removal from the wheelchair and facilitating both independent and caregiver-assisted transfers.




# Complete suspension of all bony prominences and genitals



# Tested Materials to ensure lightness, comfort and prevention

#### **Postural** inserts



#### Polyethilen foam - soft density

- · Extremely lightweight
- · Absorbs vibrations


#### Combined padding kit

Polyurethane foam of different density, lift and memory foam to ensure perfect adherence to the body and movements



#### Expanded polyurethane foam, flexible and elastic


- Amortizing
- · Supports movements



#### Antidecubitus with low sping back action poliurethane foam

- · Provides pressure distribution
- · Adapts to morphology

#### Cover



#### PU/PES combined material

- antibacterial
- · non-cytotoxic
- · non-irritant
- · resistant to virus penetration · breathable
- · antifungal
- · fire retardant

- · biocompatible
- · OEKO-TEX®
- · latex-free
- · waterproof

## **J** Senso Compacto



Shaped outer side for trochanters suspension

Removable postural inserts to adjust pelvic misalignments in the sagittal plane



The cushion is shown without the upholstery cover

Shaped inner side for ischia suspension


2 separate pelvis and femurs support elements

Empty central part for sacrum, coccyx and genitals suspension

**SENSO COMPACTO** is the **Postural Cushion** of the SENSO Range with a simple but innovative structure.

The original contoured shape of the two separate modular padded pelvis and femur support elements ensures complete relief of the areas exposed to greater pressure

The outer and inner side of each element is shaped for perfect suspension of the trochanters and ischia; The empty central part between the two elements ensures sacrum, coccyx and genitals suspension.



Size (H 7cm)	3638	3642	3840	3845	4040	4045	4242	4245	4250	4545	4550	4848	4850	5050
Width	36	36	38	38	40	40	42	42	42	45	45	48	48	50
Depth	38	42	40	45	40	45	42	45	50	45	50	48	50	50
Effective wheelchair external tubes width (indicative values)	35/38	35/38	37/40	37/40	39/42	39/42	41/44	41/44	41/44	44/47	44/47	47/50	47/50	49/52

**SENSO COMPACTO** can be customized by choosing the color of the zippers!









# SensoInpocket علا



2 separate shaped, modular, and padded trunk support elements

Empty central part for spine suspension



Length and tension adjustable strap system

**SENSO INPOCKET** is the **Postural Backrest** of the SENSO Range with a simple but innovative structure.

Two separate **modular** padded trunk support elements with a number of different-shaped postural inserts inside designed to **accommodate and/or address** the shape of the user's spine.



Size	3630	3635	3640	4030	4035	4040	4430	4435	4440	4830	4835	4840
Backrest height	30	35	40	30	35	40	30	35	40	30	35	40
Backrest width	31	31	31	36	36	36	39	39	39	43	43	43
Effective wheelchair external tubes width (indicative values)	35/38	35/38	35/38	39/42	39/42	39/42	43/46	43/46	43/46	47/50	47/50	47/50









Manufactured and distributed by:



Pro Medicare srl Via A. Montagna - Z.I. - Lotto 41 72023 Mesagne (Br) - Italy tel. (+39) 0831 77 78 40 fax (+39) 0831 73 07 39 info@promedicare.it www.promedicare.eu

- **f** ProMedicareIt
- @\_promedicare\_
- promedicare



#### 2 YEARS WARRANTY



The manufacturer reserves the right to make changes without notice. The contents of this catalogue are property of Pro Medicare S.r.l. - any improper or unauthorized use is strictly forbidden

CERTIFIED COMPANY WITH QUALITY
MANAGEMENT SYSTEM ISO 9001/ ISO 13485



